

For detecting seizure-like symptoms

Tracking symptoms typical of absence seizures and tonic-clonic seizures



Specially adapted Apple Watch

- Apple quality & value
- Robust construction
- Medical-grade sensors for heart rate, oxygen saturation, heart rhythms

Customizable for common symptoms

- High or low heart rate (sampled often)
- Dropping oxygen saturation (sampled sporadically throughout day or night)
- Convulsive motion detection*
- Trend analysis via web portal

Caregiver Apps

- Both iOS and Android
- Optional Web portal
- Location and wellness at a glance
- Highly configurable



Personalizable Watchband colors



Note that BoundaryCare is not a medical device and is not FDA cleared. It cannot diagnose or treat any medical condition.

www.boundarycare.com/order

Specifications

Dimensions	1.57 x 1.34 x 0.42 inches
Weight	.98 ounce
Connectivity	Cellular + Wifi
High heart rate alert	5 BPM increments
Low heart rate alert	5 BPM increments
A-Fib alert	Standard
Oxygen saturation alert	1% units (requires upgrade)
Convulsive motion*	Customizable
Hard fall alerts	Within 1 minute, with location
Soft fall alerts*	From 2-3 minutes
Health metric alerts	From 2-3 minutes**
Feedback to wearer	Yes, customizable
Calling to watch	With auto-answer (if desired)
PERS Emergency button	On screen (if desired)
Alerts via...	App, Text, and Email
Alert structure	Highly customizable
Night alerts	Loud & persistent
Number of caregivers	Unlimited
Caregiver Apps	iPhone, Android & Web
Water resistant	Down to 160 feet
Battery Life	24 hours or more
Protective case	Standard
Tamper-resistant band	Available (small upcharge)
Ultimate Kit	Built on Apple Watch Ultra

Additional features

Locations	GPS + cell/wifi triangulation
Location history	Last 12 hours
Watch removed alert	Standard
Low battery alert	Standard
Disconnection alert	Standard
Charging reminders	Schedulable
Customizable reminders	Unlimited, schedulable
* movements associated with tonic-clonic seizures	** refresh frequency set by user; O2 readings are occasional